

acute momentum

Trendelenburg



kwalu.

Furnishing the Future™

acute momentum

Trendelenburg

Specifically designed for acute care patient and treatment rooms, this is the Acute Momentum Trendelenburg. This chair features the quick access Trendelenburg foot pedal. Push on arms to engage reclining mechanism; footrest and back work in sync together it is not possible to independently control just the footrest.



Model: 880502
28.25W 32D 43H 22.75SW 20.5SD 19SH 25.5AH
66 Ext. Depth

Features



Nightingale Award

Kwalu's GetWell Patient Room was awarded a Nightingale Innovation Award for Surfacing Material



Durable Finish

Kwalu's solid surface (1/8" thick) proprietary finish doesn't have any of the drawbacks of wood. The frames are moisture impervious, graffiti-resistant, easy to clean and maintenance-free.



Easy to Grasp Handgrip

Extended arm creates easy to grasp handgrip.



Trendelenburg Mechanism


While at the fully reclined position depth of 66 inches only staff can engage the latch accessible at the base on either side of the chair while pushing on back to activate the Trendelenburg mechanism.



Easy to Clean

Clorox Healthcare and Kwalu are partnering to fight HAI's. Kwalu chairs have no topical coatings, so sanitizing solutions commonly used in healthcare environments can be used – including 100% bleach – without degrading the finish.

Statement of Line



880502

28.25W 32D 43H 25.5AH
66 Ext. Depth

kwalu | Furnishing the Future™

6160 Peachtree Dunwoody Rd. Building C, Atlanta, GA 30328

1-877-MY-KWALU kwalu.com

