

# modena

---

Patient



kwalu.

Furnishing the Future™

# modena


## Patient

The Modena Patient Chair with Flex Back is designed with a flex frame connecting the seat to the back. This frame enables a modest amount of flexibility in the back only, which adds a level of comfort while waiting. The seat, arms and legs remain stationary providing ergonomic support as the user changes position while unique tapered arms with a faceted edge facilitate ingress and egress. This chair works well in healthcare waiting spaces and patient rooms. It has wallsaver back legs to ensure the walls remain in pristine condition.



Model: MDAP11  
24W 27D 43H 20SW 18.5SD 19SH 26AH

## Features



**Best of NeoCon Gold**  
Kwalu's Modena Patient Chair wins gold in prestigious Best of NeoCon Awards competition.



**Nightingale Award**  
Our Modena Patient Chair was awarded a Nightingale Innovation Award for Surfacing Material.



**Faceted-Edge Frame Detailing**  
The outside face of the chair frame showcases faceted-edge detailing reminiscent of hand-crafted wood workmanship.



**Flex Back**  
Frame enables a modest amount of flexibility in the back only; seat, arms and legs remain stationary.



**Wallsaver**  
Wallsaver back legs help keep walls in pristine condition.



**Durable Finish**  
Kwalu's solid surface (1/8" thick) proprietary finish doesn't have any of the drawbacks of wood. The frames are moisture impervious, graffiti-resistant, easy to clean and maintenance-free.



**Steel Reinforced**  
Our patented steel-reinforced joint system uses over 20 pieces of steel in most Kwalu chairs. Stop ongoing and costly furniture replacement.



**Visible Cleanout**  
Cleanout is incorporated behind the seat cushion to prevent buildup of debris while providing a clean look.



**Back Cushion and Seat Cushion are Field Replaceable**  
Both the back cushion and the seat cushion can be replaced in the field.



**Brake Under Load Casters**  
Optional Brake Under Load Casters lock when a person sits on the chair, and unlock when they're not for easy movement of the chair.



**Brake Under Load Casters – Locking and Non-Locking**  
2 front locking and 2 rear non-locking weight-activated casters allow for easy movement of the chair when a person is not sitting on the chair.

**kwalu** | Furnishing the Future™

6160 Peachtree Dunwoody Rd. Building C, Atlanta, GA 30328

1-877-MY-KWALU [kwalu.com](http://kwalu.com)

